



10 minute

48 hours of glowing skin

provides skin with 48 hours of continuous hydration

facefit

time:



1 MIN.



3 MIN.



3 MIN.

directions:

- 1 Cleanse the skin by applying 2 pumps of **PreCleanse** to damp esthetic wipes.
- 2 Mix 1 pump **Special Cleansing Gel** or **prescribed cleanser** with water in a cup and apply with fan masque brush. Lightly manipulate then remove with damp esthetic wipes.

Pro Product Boost: Add 1 pump of **Exfoliant Accelerator 35** to the cleanser.

Modality Option: Use an ultrasonic brush or blade on exfoliation mode for a deep cleanse.

- 3 Dispense **Skin Hydrating Booster** and press onto skin prior to or mixed with **Skin Hydrating Masque**. Apply with a fan masque brush. Remove with damp esthetic wipes.

Modality Option: Use an ultrasonic blade on penetration mode to optimize hydration.

48 hours of glowing skin (continued)

time:

directions:



4 Spritz **Multi-Active Toner** and apply **Skin Smoothing Cream**.



5 Apply **Skin Defense Booster** or prescribed **Daylight Defense**.

1 + 2
prescription

1 - Priority Product

Skin Smoothing Cream
48 hours of hydration with
Active HydraMesh technology.

2 - Secondary Products

Skin Hydrating Booster
oil-free mix in gel boosts hydration.

Special Cleansing Gel
soap-free foaming gel cleanser
for all skin conditions.